


**Short Scalpel Handle #7**

**Standard Scalpel Handle #7**

**Scalpel Blade Remover**

For safe removal of scalpel blades, simply insert scalpel with blade into the Scalpel Blade Remover and press the sides together. This will dislodge the blade from handle and hold the blade securely for release.



6.5 cm

**No. 10002-00**

12 cm

**No. 10007-12**

16.5 cm

**No. 10007-16**

**Scalpel Blade Remover and Disposal Box**

This plastic device will remove used scalpel blades safely and retain them for safe disposal. Can be used for blades that fit the #3, #4, #7 and Gross Anatomy handles.

7.5 x 5 x 2.5 cm  
Capacity approximately: 80 blades  
Carton of 10 boxes



**No. 10002-50**

**Your fingers will never come in contact with the blade.**

