


Dumont #5SF

The finest tips ever produced by Dumont.




Super fine

11 cm
0.025 x 0.005 mm
Inox

No. 11252-00

**Dumont #5
Fine Tips**

Finished under magnification and polished on the inside, the tips close flat to 3 mm and are less likely to open under excess pressure.




Biology

11 cm
0.05 x 0.01 mm
Inox

No. 11254-20

**Dumont #5
Ceramic Coated**

These forceps are non-reflective. The scratch and corrosion resistant coating guarantees durability.




Biology

11 cm
0.05 x 0.01 mm
Ceramic coated Inox

No. 11252-50

Dumont #5CO

These forceps have been honed down to make this the lightest and most precise Dumont forceps available.



Biology

11 cm
0.05 x 0.01 mm
Dumostar

No. 11295-20

FORCEPS/TWEEZERS