

**Dumont WA #5**


These specially designed forceps offer improved stability with an ergonomic body that allows the forceps to fit more comfortably in your hand. A domed end with no soldering grooves allows the forceps to be thoroughly and easily cleaned, preventing biofilm residues.

**Dumont #5XL**

FORCEPS/TWEEZERS



Side View




**Standard**

12 cm  
0.10 x 0.06 mm

Inox

**No. 11251-21**




**Biology**

12 cm  
0.05 x 0.02 mm

Inox

**No. 11252-21**



**Standard**

15 cm  
0.1 x 0.06 mm

Inox

**No. 11253-10**