


**Dumont #L5**

Slide sleeve keeps tips together so specimens may be securely held and retained. Slide may be locked in place so that forceps can be used in regular fashion.



**Standard**


11 cm  
0.1 x 0.06 mm

Inox

**No. 11253-27**

**Dumont #5  
Teflon® Coated**

These antimagnetic forceps are ideal for use in cryogenic environments and for EM samples. Teflon® coating provides a non-stick surface.



**Standard**


11 cm  
0.1 x 0.06 mm

Teflon coated Inox

**No. 11626-11**

**Dumont #5/45**

45° tips for improved visibility and control.



**Standard**

11 cm  
0.1 x 0.06 mm

Dumoxel


**No. 11251-35**

**Dumont #5/45  
Cover Slip Forceps**

Designed to pick up cover slips out of multi-well plates. The tip width is squared for a better grip.



300%



**1 mm from tip to bend**

11.5 cm  
Upper tip: 0.1 x 0.08 mm,  
Lower tip: 0.1 x 0.02 mm  
Dumoxel

**No. 11251-33**

FORCEPS/TWEEZERS