


**Dumont #7
Fine Tips**

Finished under magnification and polished on the inside, the tips close flat to 3 mm and therefore less likely to open under excess pressure.



Biology


11.5 cm
0.07 x 0.03 mm
Inox

No. 11274-20

Dumont #7



11272-40




Standard

11.5 cm
0.17 x 0.1 mm
Dumoxel

No. 11271-30

Dumostar

No. 11297-00



Biology

11.5 cm
0.07 x 0.04 mm
Dumoxel

No. 11272-30


Dumostar

No. 11297-10

Titanium

No. 11272-40

**Dumont #7
Ceramic Coated**



Biology

11.5 cm
0.07 x 0.04 mm
Ceramic coated Inox

No. 11272-50

Dumont #7B



Standard

11.5 cm
0.17 x 0.1 mm
Inox

No. 11270-20

FORCEPS/TWEEZERS