


Dumont AA




Standard

12.5 cm
0.4 x 0.2 mm
Inox

No. 11210-20

**Dumont AA
Epoxy Coated**

Slip-proof acid and electric shock resistant epoxy coating covering 17.5 mm from the tips.




Standard

12.5 cm
0.4 x 0.2 mm
Epoxy coated Inox

No. 11210-10

Dumont SS


The tips are smooth on the inside for a gentle grip. Applying full strength will not separate tips and a strong hold on the smallest specimen is ensured.



Standard

Straight
13.5 cm
0.2 x 0.12 mm
Inox


No. 11203-23



Standard

Angled 45°
13.5 cm
0.2 x 0.12 mm
Inox

No. 11203-25



Biology

Straight
13.5 cm
0.1 x 0.06 mm
Dumoxel

No. 11200-33

FORCEPS/TWEEZERS