


Balanced #5 Forceps


These forceps feature an ergonomic design and ultra fine tips for optimal comfort and accuracy. The balanced action of the forceps ensures seamless efficiency during prolonged tasks, providing ease in rotational pivoting with the rounded handle and a counterweight technology providing enhanced control and stability.



Straight

15 cm
0.1 x 0.06 mm

No. 11851-20




Straight

15 cm
0.05 x 0.02 mm

No. 11852-20

Gold Plated Forceps


Gold-plating gives corrosion resistance to most chemicals, salts and acids. These forceps are perfect for any application where corrosion resistance and toughness are primary requirements.



Straight

11 cm
0.1 x 0.06 mm
Gold plated
Stainless Steel


No. 11915-11



Curved

12 cm
0.2 x 0.1 mm
Gold plated
Stainless Steel

No. 11917-12



Angled 30°

11 cm
0.1 x 0.06 mm
Gold plated
Stainless Steel

No. 11911-11

FORCEPS/TWEEZERS