

**Forceps with Replaceable Tips**

Stainless steel body with replaceable plastic tips. Tips are made of polyvinylidene fluoride (PVDF) and are electrostatic discharge (ESD) safe. Body and plastic tips are corrosion and abrasion resistant.




Pkg. of 2 tips and 3 screws  
Replacement plastic tips

Straight / Fine

4 cm  
0.6 x 0.6 mm

**No. 11988-01**



Pkg. of 2 tips and 3 screws  
Replacement plastic tips

Straight / Very fine

4 cm  
0.5 x 0.5 mm

**No. 11988-02**




Pkg. of 2 tips and 3 screws  
Replacement plastic tips

Curved / Very fine

4 cm  
0.5 x 0.5 mm

**No. 11988-03**




**Straight**

13 cm  
Stainless Steel body  
with plastic tips  
0.5 x 0.5 mm

**No. 11980-13**

**Fine Forceps Mirror Finish**


The last 5 mm on the inside of the tips are highly polished to help prevent tissue from sticking. Ideal for embryos, nerves or other adhesive tissue.



**Straight**

11 cm  
0.1 x 0.06 mm

**No. 11412-11**




**Straight**

11 cm  
0.05 x 0.02 mm

**No. 11413-11**

**Fine Forceps Broad Shanks**

Straight forceps with sturdy coarse tips.



**Straight**

12 cm  
0.2 x 0.12 mm

**No. 11445-12**

FORCEPS/TWEEZERS