


Moria Ultra Fine Forceps

These forceps have the most finely executed points available. They are extremely light, superbly balanced and have a gentle spring action. The shanks are specially designed with an arch which ensures maximum force to the tips. The tips close flat for 6 mm from the points to give a large holding area. The tips have been expertly finished under magnification.



Moria MC40




Straight

11 cm
0.05 x 0.01 mm

No. 11370-40

Moria MC40/B




Curved

11 cm
0.05 x 0.01 mm

No. 11370-42



Moria 9980




Straight

13 cm
0.1 x 0.06 mm

No. 11399-80

Moria 9987



Curved

13 cm
0.1 x 0.06 mm


No. 11399-87

Moria Iris Forceps

Balanced with gentle spring action and tips that are perfectly matched.



Moria MC32/B




Straight
Smooth

10 cm
0.5 x 0.3 mm

No. 11373-22

Moria MC31/B




Curved
Smooth

10 cm
0.5 x 0.3 mm

No. 11373-12


Moria MC32



Serrated

No. 11370-32

Moria MC31



Serrated

No. 11370-31