

Fine Forceps


SuperGrip Forceps

The tips are processed with MicroCrystals®, a method developed by S&T® to infuse Tungsten Carbide on stainless steel. The formation of this ultra-fine coating guarantees atraumatic handling of tissue while improving the grip dramatically. The recommended use is where forceps are utilized as needle holders or when slippery tissues might be encountered.



300%
The gripping surface is similar to a match head.


FORCEPS/TWEEZERS



Straight Standard

11 cm
0.3 x 0.25 mm


No. 00108-11



Angled 45° Standard

11 cm
0.3 x 0.25 mm


No. 00109-11



Straight SuperGrip

11 cm
0.3 x 0.25 mm

No. 00632-11



Angled 45° SuperGrip

11 cm
0.3 x 0.25 mm

No. 00649-11