



FINE FORCEPS

Bonn

The 1x2 teeth forceps No. 11084-07 are the finest of their kind with 0.1 mm wide teeth. The smooth forceps No. 11083-07 have a narrow 5 mm long platform. Both forceps have a pin inside the shanks to limit the applied pressure.




20x magnification



Straight
Smooth with platform

7 cm
Tip width: 0.3 mm
Tying platform: 5 mm

No. 11083-07




Straight
1x2 Teeth

7 cm
Tip width: 0.3 mm
Teeth width: 0.1 mm


No. 11084-07

Mirror Finish Bonn Forceps

The fine Bonn forceps features a tying platform and is also available with additional 1x2 teeth. It's designed for tissue manipulation and suturing. The wide handle ensures safe handling. The fine and precise tips make it suitable for work under a high magnification.




20x magnification



Straight
Tying platform with
1x2 teeth

10.5 cm
Tip width: 0.12 mm
Tying platform: 5 mm

No. 11090-10




Straight
Smooth with platform

10.5 cm
Tip width: 0.3 mm
Tying platform: 5 mm


No. 11093-10

Mirror Finish Suture Tying Forceps with Teeth

A unique tip style fusing a suture tying platform and a 1x2 teeth tip for tissue manipulation.



20x magnification



Straight
Tying platform with
1x2 teeth

9 cm
Tip width: 0.3 mm
Tying platform: 6 mm

No. 11090-09