


**Sliding Lock Forceps**




**Straight  
Smooth**

11 cm  
1 x 0.75 mm

**No. 11062-10**

**Episcleral Forceps**

The tips of these forceps arch toward each other, closing flat. Exceptional gripping power is possible with minimal damage to tissue. A limit pin prevents excessive applied pressure.



300%




**Straight  
Arched**

10 cm  
1.3 x 3 x 0.2 mm

**No. 11071-10**

**Pierson Fixation Forceps**

For exceptional gripping power with no damage to sclera, vascular sheath and epineurium.



1800%

Tissue retention  
hole of 0.3 mm  
behind tips.



**Straight  
0.3 mm Hole**

13 cm  
1 x 0.5 mm

**No. 18155-13**

FORCEPS/TWEEZERS