


STANDARD FORCEPS

Plastic Forceps


These forceps are non-magnetic, acid resistant and able to withstand temperatures up to 240 °C (464 °F). Tips do not spread under pressure. Plastic forceps are made of Acetyl and will be sold in packages of 10 pieces.



Straight Sharp

11.5 cm
0.4 x 0.25 mm


No. 11700-00



Straight Blunt

11.5 cm
1.5 x 0.4 mm


No. 11700-03



Curved Sharp

11.5 cm
0.3 x 0.2 mm


No. 11700-01



Angled 45° Sharp

11.5 cm
0.3 x 0.2 mm

No. 11700-02




Curved Blunt

11.5 cm
1.5 x 0.4 mm

No. 11700-04

Cross Action Forceps

These tissue transport forceps have a gentle action.




Angled Cross serrated

10 cm
4 x 8 mm

No. 11075-00


Cover Glass Forceps



Straight Smooth

10.5 cm
4 x 10 x 0.6 mm

No. 11074-02



Angled Smooth

10.5 cm
4 x 10 x 0.6 mm

No. 11073-10