

Adson DeBakey

Ideal for gentle, atraumatic gripping of vessels, tissue and other delicate structures. Two rows of atraumatic teeth interfit with one opposite row of atraumatic teeth.



Atraumatic DeBakey Serrations

This tooth pattern, as shown, consists of one row of atraumatic teeth meshed with two opposite rows of atraumatic teeth.

These two forceps are ideal for gentle, atraumatic gripping of vessels, tissue and other delicate structures.

Adson-Graefe



Adson-Brown



FORCEPS/TWEEZERS

Straight DeBakey serrations

12.5 cm
1.75 x 1.5 mm

No. 11617-12

Straight DeBakey serrations

15 cm
1.75 x 1.5 mm

No. 11616-15

Straight Fine teeth

12.5 cm
4.5 x 1.25 mm

No. 11030-12

Straight Shark teeth

12 cm
2 x 1.4 mm

No. 11627-12