

**Semken**

The Semken tweezers are characterized by their slim shape and fine tips. Their anatomical (toothing) or surgical (1x2 teeth) tips make them ideal for grasping tissue. The serrated grip surface ensures a secure hold during use.



11608-15



**Straight Serrated**

13 cm  
Tip width: 1.3 mm

**No. 11008-13**

15 cm


Tip width: 1.5 mm

**No. 11008-15**

15 cm

Tip width: 1.5 mm  
Titanium

**No. 11608-15**



**Curved Serrated**

13 cm  
Tip width: 1.3 mm

**No. 11009-13**

15 cm

Tip width: 1.5 mm

**No. 11009-15**



**Straight 1x2 Teeth**

12.5 cm  
Tip width: 1.3 mm

**No. 11024-12**

**DeBakey Round Handle Forceps**

Ideal for gentle, atraumatic gripping of vessels, tissue and other delicate structures. Two rows of atraumatic teeth interfit with one opposite row of atraumatic teeth. The round handles allow easy fingertip adjustments.



**Straight DeBakey serrations**

15 cm  
Tip width: 1.2 mm

**No. 11076-15**