


# LIGATION

## V-Hook

Aids in feeding fine suture around blood vessels, works very similar to ligation aids.



300%



11 cm

Tip Dimensions: 0.4 x 0.2 mm

No. 18052-12

## Ligation Aids

A handy instrument for leading a fine suture under and around a blood vessel for tying-off purposes. Also useful for ligating peripheral nerves. Grooves toward hole help threading.



18063-13



300%



10.5 cm

Hole Ø: 0.4 mm

No. 18062-12



200%


200%



13.5 cm

Hole Ø: 1 mm

No. 18062-13



13.5 cm

Hole Ø: 1 mm

No. 18063-13

VASCULAR INSTRUMENTS