


**Mini Hippocampal Dissecting Tool**

The instrument ends are malleable to allow you to adjust their shape to the correct angle for optimal positioning.




14.5 cm  
20 x 6 mm & 13 x 3 mm  
Tip thickness: 0.5 mm

**No. 10099-12**

**Hippocampal Dissecting Tool**

The instrument ends are malleable to allow you to adjust their shape to the correct angle for optimal positioning.



16 cm  
27 x 14 mm & 27 x 6 mm  
Tip thickness: 0.5 mm

**No. 10099-15**