

Spoon

Mini Perforated Spoon

Perforated Spoon

Mirror Finish Cutting Spoon

Most commonly used to take embryos out of solution and handling very small organs.

Most commonly used to take embryos out of solution and handling very small organs.

Non stick cutting spoon.

Moria

Moria



Moria




Moria 1121B


Moria MC 17C

Moria MC 17 BIS


Moria MC 17




12.5 cm
Tip Ø: 8 mm
Depth: 2 mm



10 cm
Tip Ø: 8 mm
Depth: 1 mm
Hole Ø: 0.8 mm



14.5 cm
Tip Ø: 15 mm
Depth: 3 mm
Hole Ø: 0.9 mm



14.5 cm
Tip Ø: 20 mm
Depth: 5 mm
Hole Ø: 0.9 mm



13 cm
Tip Ø: 4.5 mm
Depth: 0.15 mm

No. 10321-08

No. 10370-19

No. 10370-18

No. 10370-17

No. 10360-13

SPOONS