

NEEDLE HOLDERS

Castroviejo

Micro pattern, with slightly blunted tips.



Never use a micro needle holder with needles that are too thick, as this can easily damage the fine jaws. Use only micro needles of less than 0.2 mm diameter (see page 154).



Straight

9 cm
Tip width: 0.6 mm
Lock: Yes

No. 12060-01

Tip width: 0.6 mm
Lock: No

No. 12060-02



Curved

9 cm
Tip width: 0.6 mm
Lock: Yes

No. 12061-01

Tip width: 0.6 mm
Lock: No

No. 12061-02

Round-Handled Micro/Macro Needle Holder

These Round-Handled Needle Holders have a dual-purpose and can be used for microscopic and macroscopic applications. The fine tapered tips can be used to handle small needles up to 0.2 mm in diameter. While the remaining portion of the tip ensures a secure hold for larger needles. The balanced design improves control and stability, while the ergonomic handle enhances grip and reduces hand strain.



Curved

13.5 cm

Lock: Yes

No. 12875-13