


**AutoClip® System**

The AutoClip® System is a safe, secure method of wound closure. Faster and easier than suture, this system is perfect when staples may be too thin to provide adequate closure. **Non-sterile. Components sold separately or as a kit.** We recommend using wet steam sterilization at 120-125 °C for at least 30 minutes in an autoclave to sterilize the device.



**Applier**



**Clips**



**Remover**

WOUND CLOSURE

*AutoClip® Kit*

- 1 Applier
- 9 mm clips (100)
- 1 Remover
- Storage case

*AutoClip® applier*

11.5 cm

**No. 12020-00**

*AutoClip® clips*

9 mm  
Pkg. of 100 clips

**No. 12022-09**

*AutoClip® clips*

9 mm  
Pkg. of 1,000 clips

**No. 12022-19**

*AutoClip® remover*

12.5 cm

**No. 12023-00**