

Marking Pen

The Marking Pen was developed for easy, long lasting identification of laboratory animals. The markers contain no harmful chemicals and instead nontoxic pigment approved by the Japanese Ministry of Health for use in food, drugs and cosmetics. It is intended for use on the fur and not the skin and will last for 6-12 weeks.



CAUTION: Color can stain skin and clothing. The use of gloves and lab coats is recommended. These colors will not work on black fur.



Black	Blue	Green	Red	Yellow	Orange
Fine Tip 4.5 x 7.5 mm 10 ml	Fine Tip 4.5 x 7.5 mm 10 ml	Fine Tip 4.5 x 7.5 mm 10 ml	Fine Tip 4.5 x 7.5 mm 10 ml	Fine Tip 4.5 x 7.5 mm 10 ml	Fine Tip 4.5 x 7.5 mm 10 ml
No. 24202-06	No. 24202-02	No. 24202-01	No. 24202-04	No. 24202-03	No. 24202-05



Black	Blue	Green	Red	Yellow	Orange
Standard Tip 15 x 10 x 10 mm 15 ml	Standard Tip 15 x 10 x 10 mm 15 ml	Standard Tip 15 x 10 x 10 mm 15 ml	Standard Tip 15 x 10 x 10 mm 15 ml	Standard Tip 15 x 10 x 10 mm 15 ml	Standard Tip 15 x 10 x 10 mm 15 ml
No. 24202-16	No. 24202-12	No. 24202-11	No. 24202-14	No. 24202-13	No. 24202-15

Surgical Skin Marker

This marker features a versatile tip for precise fine (0.2 mm) line marking, or when held at an angle, a broad (2 mm) line marking. The marker's violet ink is nontoxic, non-tattooing and surgical prep-resistant. Markings can be removed from the skin with isopropyl alcohol. A six-inch ruler calibrated in both centimeters and inches is included.



No. 18000-30

50% of actual size

ANIMAL ACCESSORIES