

Instrument Marking Tape

Organize instruments quickly and easily with this economical, non-toxic, pressure sensitive tape. Color is unaffected by boiling, cold sterilization or autoclaving. 8 colors available. **Autoclavable.**




| Color | Roll length (cm) | Tape width (cm) | Material | Item |
|---------------|------------------|-----------------|----------|---------------------|
| Green | 290 | 0.6 | Vinyl | No. 29071-01 |
| Blue | 290 | 0.6 | Vinyl | No. 29071-02 |
| Yellow | 290 | 0.6 | Vinyl | No. 29071-03 |
| Red | 290 | 0.6 | Vinyl | No. 29071-04 |
| White | 290 | 0.6 | Vinyl | No. 29071-05 |
| Black | 290 | 0.6 | Vinyl | No. 29071-06 |
| Orange | 290 | 0.6 | Vinyl | No. 29071-08 |
| Purple | 290 | 0.6 | Vinyl | No. 29071-09 |

**Instrument Marking Tape
7 assorted colors**

Organize instruments quickly and easily with this economical, non-toxic, pressure sensitive tape. Color is unaffected by boiling, cold sterilization or autoclaving. Includes 7 assorted color rolls: Yellow, Red, Black, White, Green, Blue and Orange. **Autoclavable.**



| | | | | |
|--|-----|-----|-------|---------------------|
|  | 107 | 0.6 | Vinyl | No. 29080-02 |
|--|-----|-----|-------|---------------------|

Instrument Marking Tape

Organize instruments quickly and easily with this economical, non-toxic, pressure sensitive tape. Color is unaffected by boiling, cold sterilization or autoclaving. 9 colors available. **Autoclavable.**



| Color | Roll length (cm) | Tape width (cm) | Material | Item |
|---------------|------------------|-----------------|----------|---------------------|
| Green | 510 | 0.6 | Vinyl | No. 29070-01 |
| Blue | 510 | 0.6 | Vinyl | No. 29070-02 |
| Yellow | 510 | 0.6 | Vinyl | No. 29070-03 |
| Red | 510 | 0.6 | Vinyl | No. 29070-04 |
| White | 510 | 0.6 | Vinyl | No. 29070-05 |
| Black | 510 | 0.6 | Vinyl | No. 29070-06 |
| Brown | 510 | 0.6 | Vinyl | No. 29070-07 |
| Orange | 510 | 0.6 | Vinyl | No. 29070-08 |
| Purple | 510 | 0.6 | Vinyl | No. 29070-09 |