

Colored Aluminum Cases

This instrument case comes with a spiked silicone mat for secure storage of surgical instruments. Five colors available but may vary due to lot to lot variations.

**The case and silicone insert are autoclavable but should not be used for sterilization of surgical instruments.**



Blue  
18 x 9 x 3 cm

**No. 20330-00**



Red  
18 x 9 x 3 cm

**No. 20330-02**



Green  
18 x 9 x 3 cm

**No. 20330-01**



Yellow  
18 x 9 x 3 cm

**No. 20330-03**



Silver  
18 x 9 x 3 cm

**No. 20330-04**



INSTRUMENT CARE  
& ACCESSORIES