

Utility Instrument Case

Comes with a silicone mat for safe placement of instruments. **Case and silicone insert are autoclavable but should not be used for sterilizing surgical instruments.**

Moria



Silver
18 x 9 x 3 cm

No. 20311-21

Instrument Case

This plastic case is compact, safe, secure, and slotted for placement of individual surgical instruments, aiding in the prevention of tip damage. Seven colors available. **Case is autoclavable but should not be used for sterilizing surgical instruments.**



Instruments not included.

Blue

Green


Clear

Yellow

Pink

Orange

Purple



17 x 8 x 2 cm

17 x 8 x 2 cm

17 x 8 x 2 cm

17 x 8 x 2 cm

17 x 8 x 2 cm

17 x 8 x 2 cm

17 x 8 x 2 cm

No. 20830-00

No. 20830-01

No. 20830-02

No. 20830-03

No. 20830-04

No. 20830-05

No. 20830-06