

**High Temperature Cautery Kit**

This high temperature cauterizer is a great alternative to disposable cautery units. Temperature at tips is approximately 1000 °C (1832 °F). Replacement cautery tips are available. **Not autoclavable.**



Plastic body  
 Stainless steel tips: 0.2 mm Ø  
 Operation temperature: 1000 °C (1832 °F)  
 Total weight: 70 g  
 Length: 16 cm (without tip)

**Cautery Kit**  
**No. 18010-00**

- Cautery kit includes:
- 1 Cautery handle
  - 2 AA batteries
  - 2 Standard tips
  - 2 Fine tips
  - 1 Storage case



Replacement tips  
**Single use**  
 Straight / Standard  
 Sterile: Yes

Box of 10 tips  
**No. 18010-01**



Replacement tips  
**Single use**  
 Straight / Fine  
 Sterile: Yes

Box of 10 tips  
**No. 18010-02**

**Low Temperature Cautery Kit**

This low temperature cautery kit is a great alternative to disposable cautery units. Temperature at the tips is approximately 649 °C (1200 °F). Replacement cautery tips are available. **Not autoclavable.**



Plastic body  
 Stainless steel tips: 0.2 mm Ø  
 Operation temperature: 649 °C (1200 °F)  
 Total weight: 60 g  
 Length: 12 cm (without tip)

**Cautery Kit**  
**No. 18019-00**

- Cautery kit includes:
- 1 Cautery handle
  - 1 AA battery
  - 2 Fine tips
  - 2 Elongated fine tips
  - 1 Storage case



Replacement tips  
**Single use**  
 Straight / Fine  
 Sterile: Yes

Box of 10 tips  
**No. 18019-01**



Replacement tips  
**Single use**  
 Elongated Straight / Fine  
 Sterile: Yes

Box of 10 tips  
**No. 18019-02**