

SCISSORS

Vannas Spring Scissors

Ideal for introductory courses in microdissection. The manufacturer, material and general pattern are the same as our No. 15003-08 but with much of the costly hand finishing omitted.



**Straight
Sharp/Sharp**

9 cm
Cutting edge: 5 mm
Tip Ø: 0.35 mm

No. 91500-09



**Curved
Sharp/Sharp**

9 cm
Cutting edge: 5 mm
Tip Ø: 0.35 mm

No. 91501-09

Fine Scissors



**Straight
Sharp/Sharp**

11.5 cm
Cutting edge: 16 mm

No. 91460-11



**Curved
Sharp/Sharp**

11.5 cm
Cutting edge: 16 mm

No. 91461-11

Sanvenero



**Straight
Sharp/Sharp**

14.5 cm
Cutting edge: 22 mm

No. 91408-14