


Micro Adson




80% of actual size



Serrated

12 cm
Tip width: 0.7 mm

No. 91118-12



1x2 Teeth

12 cm
Tip width: 1.3 mm

No. 91119-12

Adson

Wide, strong grip with narrow tips.




80% of actual size



Serrated

12 cm
Tip width: 1.5 mm

No. 91106-12




1x2 Teeth

12 cm
Tip width: 1.5 mm

No. 91127-12

Adson-Brown



Shark teeth

12 cm
Tip width: 2.5 mm

No. 91162-12

Adson with Cross Serrations



**1x2 Teeth/
Cross-serrations**

12 cm
Tip width: 1.3 mm

No. 91427-12